시차주사열량계 (DSC, Differential Scanning Calorimeter)

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- Glass transition temperature (유리전이온도, T_g)
- Specific heat capacity (비열)
- Temperature of fusion (용융점, T_m)
- Heat of fusion (용융열, ΔH)
- Crystallinity, degree of crystallization (결정화도)
- Melting behavior, fraction melted (용융거동)
- Temperature of crystallization (결정화 온도, T_c)
- Heat of crystallization (결정화열, ΔH)
- Cold crystallization (재결정화)
- Chemical reaction (화학반응)
- Curing (경화), degree of cure (경화도), 경화열 (ΔH), post-curing
- Polymorphism (change of crystal modification), mesophase transition, solid-solid transition
- Evaporation, desorption (moisture), vaporization, sublimation
- Thermal decomposition (pyrolysis, depolymerization)
- Thermal stability (열안정성)
- Lambda transition
- Oxidative degradation, oxidation stability (산화안정성), Oxidation Induction Time (OIT, 산화도입기)
- Purity (순도)등의

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Part II.

2. DSC 측정 곡선 해석 (Curve interpretation)
 - 2.1. Dynamic curve interpretation
 - 2.2. Isothermal curve interpretation
3. DSC 응용 사례 (Application)
 - 3.1. Analysis of Polymers
 - 3.2. Analysis of Elastomers
 - 3.3. Analysis of Pharmaceuticals
 - 3.4. Analysis of Foods
 - 3.5. Analysis of Composites

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

Measures the difference of heat flow to sample and to reference side using one single furnace. Measures sample and reference temperature using single or multiple thermocouples (thermopile).

Measures the difference of heat flow between sample side and reference side using two separate small furnaces. Measures sample and reference temperature with platinum resistance thermometer

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- Glass transition temperature (유리전이온도, T_g)
- Specific heat capacity (비열)
- Temperature of fusion (용융점, T_m)
- Heat of fusion (용융열, ΔH)
- Crystallinity, degree of crystallization (결정화도)
- Melting behavior, fraction melted (용융거동)
- Temperature of crystallization (결정화 온도, T_c)
- Heat of crystallization (결정화열, ΔH)
- Cold crystallization (재결정화)
- Chemical reaction (화학반응)
- Curing (경화), degree of cure (경화도), 경화열 (ΔH), post-curing
- Polymorphism (change of crystal modification), mesophase transition, solid-solid transition
- Evaporation, desorption (moisture), vaporization, sublimation
- Thermal decomposition (pyrolysis, depolymerization)
- Thermal stability (열안정성)
- Lambda transition
- Oxidative degradation, oxidation stability (산화안정성), Oxidation Induction Time (OIT, 산화도입기)
- Purity (순도) 등의 열적 거동 (Thermal effects)을 분석할 수 있다.

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Copyright© YEONJIN Corp. Scientifics